

# Formfit® Pro Knee

Runner's Knee?

Pain when squatting?

Concern or Instability?

Össur Formfit® Pro Knee is a compressive knee sleeve that features a SmartPad™ buttress for targeted support and cushioning.

The patented design realigns the patella (kneecap), which can be beneficial if you're suffering from patellar tendinitis (Jumper's Knee) or chondromalacia patellae (Runner's Knee).

**Is the Formfit® Pro Knee for you?**

**For knee conditions that may benefit from compression, such as;**

- Patella tendinitis
- Chondromalacia Patella
- Mild sprains and strains
- Knee pain
- Feeling of instability



Visit [jointpainfreedom.co.uk](http://jointpainfreedom.co.uk) to learn more or scan the QR code.

**Support and pain management,  
giving you confidence during activity**

- MotionTech™ for compression and dynamic fit with superior moisture wicking and breathability, keeping your knee cool during wear
- Certified medical grade compression (20-30 mmHg)
- Lightweight CoolVent™ knit in popliteal area for added breathability and comfort
- Patented SmartPad™ for targeted support and cushioning
- Patella graphics for intuitive placement
- Patented FoldLock™ system - easier and secure application
- Lightweight, removable spiral steel stays





# FORMFIT PROKNEE

## DESIGNED TO MAINTAIN ACTIVITY

Can be complimentary to knee management strategies that also include injection therapy

**If your knee is inflamed or painful and you'd like to understand more about fast and long-term pain relief:**

→ **Speak to your clinician about Monovisc or Cingal**

Bracing and injection therapy often work best together. The brace helps stabilise and protect your knee, while the hyaluronic acid injection supports the joint from within. Compression and off-loading can enhance the benefit of the injection and give you confidence to return to daily activity.

The injection then continues to nourish and lubricate the joint for longer-term comfort and joint health.

**89% of patients with mild to moderate arthritis have significant pain reduction after 1 week with Cingal.**



### Hybrid Cingal & Hyaluronic Acid Matrix Monovisc



- **Cingal** combines fast acting steroid with **Monovisc**, a proven high dose, high performance hyaluronic acid.
- Studies taken from 3 separate **Cingal** level 1 human trials show:

Significant  
Pain relief  
at 1 week

Pain relief  
that lasts  
beyond 6  
months

More  
effective &  
safer vs  
steroid alone

92/100  
patients  
respond +ve  
to treatment



- **Monovisc** as a standalone treatment is a high performing, long lasting option without steroid.
- A hyaluronic acid with a range of molecular weights, that offer **Lubrication, Cushioning, Synovial Health**